

EDAN



C3A/C6A

Video Colposcope



C3A/C6A

Video Colposcope



EDAN Edan Instruments, Inc.

■ 3/F - B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
■ Tel: +86-755-26898326 ■ Fax: +86-755-26898330 ■ www.edan.com.cn ■ Email: info@edan.com.cn

■ All rights reserved. This catalogue is only a temporary version. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



ENG-OB & GYN-C3A&C6A-V1.1-20140415

Care for Health



With its bright LED cold lighting system, high-resolution camera and electronic green filter, C3A/C6A Video Colposcope realizes full-screen display of precise cervical examination, offering crucial and accurate information for clinical decision making to the gynecologist.

Features

- LED cold lighting system reproduces an original view of tissue
- Advanced and fast auto-focusing system ensures clear image
- Remote capture control and user-friendly keyboard layout for easy operation
- Electronic green filter enhances vascular image display without compromise of illumination
- Acetic acid and iodine reaction test timer to capture the right and precise image
- Real-time display of magnification
- Powerful image display and data management software

R-way™

The R-way function observes the variation of "color" on cervix to interpret the performance of cervical diseases

Saline application – 'R for Red'
To observe the area with rich blood circulation, this reflects in red.

Acetic acid 5% – 'W for White'
To observe increase in nucleus amount, this reflects in white.

CIN2+

Green filter – 'A for Atypical'
To better observe atypical vascular, this defines as A.

Lugol's iodine solution – 'Y for Yellow'
To observe decrease in cell plasma, this reflects in yellow.



C3A/C6A Video Colposcope

Patient management software

