

Eliminación el exceso de drogas y venenos



HA230

Cartucho de Hemoperfusión Desechable

18 años de evidencia clínica
Ampliamente utilizado en 60 países

HA230 Cartucho de Hemoperfusión Desechable Propiedades

Volumen de adsorbente (ml)	230±3
Volumen de sangre (ml)	145±5
Material adsorbente	Copolímeros de estireno-divinilbenceno reticulados secundarios
Material de la carcasa	Polycarbonato
Método de esterilización	Esterilización por irradiación
Paquete unitario	236 mm (largo) 100 mm (ancho) x 102 mm (alto) 0,7 kg

Jafron: el mayor fabricante de columnas de adsorción del mundo




Jafron I+D y centro de fabricación



CE



ISO 9001



ISO 13485

Referencias

[1] Jiang Tianlan. Chin J Ethnomedicine and Ethnopharmacy 2012;130-131.

[2] Yu Dan. J Chin foreign medical treatment 2012; No.02:67

[3] Xiangming Yang, et al. American J Emergency Medicine 2017

[4] Deven Juneja, et al. Indian J Critical Care Medicine 2011

[5] Shuangxin, et al. J Diabetes Investig 2018; 9:975-978

[6] Xiaoyun Si, et al. PLOSONE 2015; 10(7)

[7] Guo Zhengke, et al. Chin J Blood Purification 2002; Vol.1,No.6: 15

[8] Weng Guirong, et al. Chin Crit Care Med 2012; Vol. 24, No. 6:368

[9] Jiling Huang, er al. Medicine 2017;96:27

[10] Hui Dong, et al. Medicine 2017; 96:25

[11] Yuning Shi, et al. PLOSONE 2012; 7(7)

[12] Chen Zhi, et al. Chin J Emerg Med 2005; Vol.14, No.2: 140-143

[13] Liu Yadong, et al. Hebei Medicine 2012; Vol.18, No.12: 1726-1728

[14] Tuo Yanhong, et al. J Yangtze University 2009; Vol.6, No.2: 33-34

[15] Xue-hong ZHAO, et al. J Zhejiang Univ-Sci 2015; 16(8): 720-726

[16] An Li, et al. Blood Purification 2016; 42: 92-99

[17] Hui Dong, et al. Medicine 2017; 96:25