

Toma el control de la tormenta de citoquinas



HA330

Cartucho de Hemoperfusión Desechable

CE 0197

18 años de evidencia clínica
Ampliamente utilizado en 50 países
Sobre 30000 tratamientos anuales

HA330 Propiedades

Volumen de adsorbente (ml)	330±3
Volumen de sangre (ml)	185±5
Material adsorbente	Copolímeros de estireno-divinilbenceno reticulados secundarios
Material de la carcasa	Polycarbonato
Método de esterilización	Esterilización por irradiación
Paquete unitario	290 mm (largo) x 105 mm (ancho) x 105 mm (alto) 0,7 kg


Jafron: el mayor fabricante de columnas de adsorción del mundo




Jafron I+D y centro de fabricación



CE



ISO 9001



ISO 13485

Referencias

1.Ankawi, G. et al. A New Series of Sorbent Devices for Multiple Clinical Purposes: Current Evidence and Future Directions. Blood Purif. 47, 94–100 (2019).

2.Pomarè Montin, D. et al. Biocompatibility and Cytotoxic Evaluation of New Sorbent Cartridges for Blood Hemoperfusion. Blood Purif. 46, 187–195 (2018).

3.Huang, Z., Wang, S. R., Su, W. & Liu, J. Y. Removal of humoral mediators and the effect on the survival of septic patients by hemoperfusion with neutral microporous resin column. Ther. Apher. Dial. 14, 596–602 (2010).

4.Huang, Z., Wang, S. R., Yang, Z. L. & Liu, J. Y. Effect on extrapulmonary sepsis-induced acute lung injury by hemoperfusion with neutral microporous resin column. Ther. Apher. Dial. 17, 454–461 (2013).

5.Sun, S. et al. High-volume hemofiltration plus hemoperfusion for hyperlipidemic severe acute pancreatitis: A controlled pilot study. Ann. Saudi Med. 35, 352–358 (2015).

6.Wang, Y. T. et al. Effects of hemodialysis and hemoperfusion on inflammatory factors and nuclear transcription factors in peripheral blood cell of multiple organ dysfunction syndrome. Eur. Rev. Med. Pharmacol. Sci. 20, 745–750 (2016).

7.Li, M. qin et al. Hemodiafiltration Combined with Resin-Mediated Absorption as a Therapy for Hyperlipidemic Acute Pancreatitis. Cell Biochem. Biophys. 69, 699–702 (2014).

8.Arslan, B., Kucukbingoz, C., Kutuk, M. & Gunduz, H. A single-center experience with resin adsorption hemoperfusion combined with continuous veno-venous hemofiltration for septic shock patients. Med. Sci. | Int. Med. J. 1 (2018). doi:10.5455/medscience.2018.07.8950

9.Chavez, J. R. et al. A case of leptospirosis with acute respiratory failure and acute kidney injury treated with simultaneous extracorporeal membrane oxygenation and haemoperfusion. BMJ Case Rep. 12, 1–6 (2019).

10.Hui, L., Zeng, A., Zhang, X., Zang, K. & Shang, F. Evaluation of the therapeutic effect of hemopurification in hyperlipidemic severe acute pancreatitis. 12, 1004–1010 (2019).

11.Bai, M. et al. Continuous venovenous hemofiltration combined with hemoperfusion for toxic epidermal necrolysis: a retrospective cohort study. J. Dermatolog. Treat. 28, 353–359 (2017).

12.Yuan, H., Chen, S., Hu, F. & Zhang, Q. Efficacy of Two Combinations of Blood Purification Techniques for the Treatment of Multiple Organ Failure Induced by Wasp Stings. Blood Purif. 42, 49–55 (2016).

13.Li, Z. et al. Effects of hemodialysis combined with hemoperfusion on severe acute pancreatitis. Turkish J. Gastroenterol. 29, 198–202 (2018).

14.Tang, Y., Zhang, L., Fu, P., Kang, Y. & Liu, F. Hemoperfusion plus continuous veno-venous hemofiltration in a pregnant woman with severe acute pancreatitis: A case report. Int. Urol. Nephrol. 44, 987–990 (2012).

15.Liang, M. J. & Zhang, Y. Clinical analysis of penicillidene hydrochloride combined with hemoperfusion in the treatment of acute severe organophosphorus pesticide poisoning. Genet. Mol. Res. 14, 4914–4919 (2015).

16.Xu, X. et al. Effect of HA330 resin-directed hemoabsorption on a porcine acute respiratory distress syndrome model. Ann. Intensive Care 7, (2017).



JAFRON BIOMEDICAL CO., LTD.

Address: No.98 Technology Six Road, High-tech Zone, Zhuhai City, 519085 Guangdong, China
Tel: 0086-20-88527536; 0086-756-3619196
Email: overseatrade@jafron.com
Website: en.jafron.com